


A2200 Release Film

- A2200 is a low cost, very flexible release film which is ideal for resin infusion processing.
- Coloured opaque orange to allow resin front to be monitored during infusion plus aiding complete removal whilst de bagging.
- Suitable for use with epoxy, polyester and vinylester resins although we strongly advise that a test piece is produced prior to committing to a full scale infusion in order to determine compatibility and performance meet desired standards.
- Two standard perforation patterns available (others available on request):-

P16: Very fast resin transfer. Ideal for placing underneath resin feeder lines to enable their removal from a peel plied surface. Due the large number of perforations, the film may become a little fragile following infusion and difficult to remove from the laminate surface in certain applications and we recommend a test piece be produced before applying to a large surface area. Nominal hole diameter 0.5mm.

P31: Fast resin transfer. A standard perforation style for resin infusion. Nominal hole diameter 0.8mm.



PHYSICAL PROPERTIES

- Maximum Use Temperature
- Perforations

120°C
P16, P31

AVAILABILITY AND PACKAGING

- Thickness
- Standard Widths (others available to order)
- Standard Roll Length

0,025mm
100mm & 1300mm
500 lin m

STORAGE AND HANDLING

- Store in original packaging
- No handling problems experienced

Aerovac Systems Ltd.

Tel: 01274 550500 Fax: 01274 550501 e-mail: sales@aerovac.com website: www.aerovac.com

All statements, technical information and recommendations contained in this publication are based on tests believed to be reliable, but their accuracy and/or completeness are not guaranteed. The user shall determine the suitability for this particular purpose and shall assume all risk and liability in connection herewith.

revision: N/A new